	<p>Gebr. Brasseler GmbH &amp; Co. KG Trophagener Weg 25 32657 Lemgo - Germany</p> <p>Telefon +49 (0) 526 1 701-400 Telefax +49 (0) 526 1 701-580 info@kometmedical.de www.kometmedical.de</p>
	0197

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**(EN)**

Instructions for use and safety recommendations for the application of rotary and reciprocating instruments in the medical field

The instruments are made of high-grade materials. They are designed, constructed and produced with utmost care. These quality instruments assure best working results and long service life provided they are used in the proper manner. Therefore, the following instructions for use and safety recommendations have to be observed. Improper use can lead to damage to the tissue, premature wear, destruction of the instruments and injury to the operator, patient or other persons.

1. Proper use

- Observe this information and the instructions for use for surgical motor systems when using rotary and reciprocating instruments.

- Prior to use, make sure that the instruments are in perfect technical condition and that sterility is assured.

Damaged instruments should not be reused.

- Insert all rotary instruments in the hand piece until the stop position is reached.

- Instruments with clamping marks (line or ring) have to be lowered into the hand-piece until the marks are no longer visible.

- Prior to operation, make sure that the instrument is firmly positioned.

- Do not use the instruments for work on metal (e.g. implants) and do not allow them to come in contact with metal during use as this might cause nicks in the blades.

- Craniotomes are to be used with corresponding dura guard only.

- To ensure traceability of the MD by means of the UDI, please make sure to keep the packaging of articles without direct marking.

2. Indication

The instruments have been designed for working on bone, skin, soft tissue and cements. Instruments of different materials adapted to each application are available.

- Dermabrasion:

instruments with diamond-coated working part

- Maxillo-facial surgery:

instruments made of stainless steel, tungsten carbide instruments, diamond instruments

- ENT surgery:

instruments made of stainless steel, tungsten carbide instruments, diamond instruments

- Neurosurgery:

stainless steel, carbide and diamond instruments

- Orthopedics:

stainless steel, carbide and diamond instruments

The surgeon is responsible for the selection of the appropriate instrument.

3. Contraindication

The instruments are not suited for use in the CNS (brain and spinal cord).

4. Contact pressure

Excessive contact pressure has to be avoided.

- It can lead to damage of the instruments' working part (steel and tungsten carbide instruments) as well as blade breakouts. At the same time extreme heat is generated which increases the risk of thermal necrosis.

- Excessive contact pressure can lead to the abrasive grit of diamond instruments breaking off.Increased heat development will be the result.

- At worst, breakage of the instrument can occur if excessive contact pressure is applied. Jamming or leverage of the instrument must be avoided during operation.

Optimal material reduction is assured at a contact pressure of approx. 2N. This can be compared with the pressure which is applied when writing with a fountain pen.

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5. Cooling

- When using rotary and reciprocating instruments sufficient cooling has to be supplied, otherwise the instrument blades will be clogged with chips which leads to increased heat generation and at worst to irreversible damage to the bone (thermal necrosis). Moreover, the service life of the instruments will be reduced.

6. Sorting out blunt and damaged instruments

Examine instruments for bluntness and damage with magnifying glasses. Look for the following indications:

- Blunt (rounded) blades and breakouts
- Blank spots on diamond instruments
- Damage to the shank
- Bent or "out of true" instruments

Under no circumstances use blunt and damaged instruments any more. Otherwise instrument breakage and risk of injury to the patient and operator may be the result.

7. Cleaning, disinfection, sterilization and storage

- A validated reprocessing method is available on our web site www.kometmedical.de under "Instrument reprocessing".

- Prior to first use and immediately after each use, rotary instruments have to be cleaned, disinfected and sterilized. The instruments should be stored in their original packages at room temperature, protected against dust and humidity.

- The instruments should be placed in hygienic bur blocks, dishes or other suitable containers. The same applies to sterilized instruments and instruments in sterile packages. The instruments are to be stored in an environment free of dust and humidity, to exclude recontamination.

- After surgical operation, immediately remove blood, secretion, tissue and bone residues. Dried residues can cause corrosion.

Make sure that all residues are removed; this applies especially to hollow burs.

- The instruments can be cleaned in the ultrasonic bath. However, care should be taken that the instruments do not come in contact with each other; otherwise the blades can be damaged.

- It is recommended to use anticorrosive disinfection and cleaning agents (e.g. Komet DC1).

- Contact with H<sub>2</sub>O<sub>2</sub> (hydrogen peroxide), chlorous disinfection and cleaning agents as well as oxalic acids are to be avoided because these agents may lead to pitting and corrosion.

- For immersion time and suitability of disinfection and cleaning substances for specific instrument types see instructions for use of the manufacturer of these agents.

- Disinfection and cleaning agents must be rinsed off thoroughly under running water immediately after use; otherwise corrosion can occur due to chemical reactions.

- Disinfection of the instruments can be carried out in the gentle cycle (e.g. Vario TD) of the washer/disinfector. Use cleaning and disinfection agents with approved material compatibility. We recommend deconex 28 AlkaOne (co. Borer Chemie).

- Sterilization should take place in the steam sterilizer.

- A temperature above 180°C has to be avoided. Exceeding this temperature would result in loss of hardness of the working part and thus reduction of service life.

Make sure that each selected disinfection and sterilization method is suitable for the specific instruments. The instruments may only be returned in clean, disinfected and sterilized condition.

8 . Safety and liability

- The user is responsible to check the product prior to use whether it is suitable for the intended purpose.

- The user is responsible for the application of the instruments.

- In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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### EN

5. Colding
When using rotating and reciprocating instruments, sufficient cooling should be provided. Otherwise, the instrument teeth will be blocked by chip accumulation, which can lead to excessive heat transfer and will lead to irreversible bone damage (thermal death). In addition, the service life of the instrument will be shortened.

6. Treatment of worn and damaged instruments
Inspect worn instruments with a magnifying glass. Pay attention to:
- chipped or uneven teeth (brilliant or rounded)
- missing particles on diamond instruments
- damage to the handle
- bent or deformed instruments
Regardless of the situation, worn instruments must not be used again. Otherwise, the instrument will be damaged and lead to patient and operator injury.

7. Cleaning, disinfection, sterilization and storage
- A validated preparation method can be downloaded from our website: www.kometmedical.de Select product information download cleaning and disinfection method.
- Before using for the first time, clean, disinfect and sterilize the instruments.
- Until used for the first time, store the instruments in their original packaging at room temperature, protected from dust and humidity.
- The preparation of the instruments should be carried out in trays, burs or other types of containers in good hygienic conditions. The same applies to sterile instruments and instruments in sterile packaging. The storage of instruments must be carried out in conditions free of dust, humidity and recontamination.
- After the operation, the instruments must be removed immediately from the patient. Do not allow the instruments to dry! Dried residues cause corrosion. This instruction must be observed especially in tubular instruments!
- Instruments can be cleaned using ultrasonic cleaning. However, care must be taken to ensure that the instruments are not damaged.
- Only use liquid disinfectants and cleaning solutions with corrosion protection (e.g. Komet DC1).
- Avoid contact with H2O2 (hydrogen peroxide), disinfectants and cleaning solutions containing chlorine or oxalic acid. If these instructions are not followed, corrosion and/or holes can occur in the instruments.
- The form of use, immersion time and suitability of the cleaning and disinfection liquids for each instrument type depend on the manufacturer's instructions.
- Disinfection and cleaning liquids must be thoroughly rinsed off after use with running water. If not followed, chemical reactions can occur which lead to corrosion.
- Instruments must be stored in a clean, dry and dust-free environment.
- Do not use temperatures above 180°C. Overheating can lead to a loss of strength in the active part of the instruments, as well as a reduction in service life.
During disinfection and sterilization, special attention must be paid to ensure that the selected procedure is suitable for the instrument in question. Instruments must only be returned in a clean, disinfected and sterilized state.

8. Safety and responsibility
- The operator is obliged and is the only person responsible for checking - before use - that the product is suitable for the intended use.
- The use of the instruments is the responsibility of the operator.
- Negligence on the part of the operator will result, in the event of damage caused by the operator, in a reduction or complete exclusion of liability on the part of Gebr. Brasseler. This will be especially the case if our recommendations for use or the warnings were not observed, or if use was unintentionally incorrect.

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© Gebr. Brasseler GmbH & Co. KG