

Καταμετρήσεις	Μακρ. κίρροσ/κν	Κίρροσ/κν	Πίτνυμκν	Κίρροσ/κν	Κίρροσ/κν
NSXXXX.025	20.000 rpm	Aesculap Low Speed	ei olenn	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 rpm	Aesculap Low Speed	ei olenn	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 rpm	Aesculap Hilan	1 rengas	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 rpm	Aesculap Hilan	2 rengas	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 rpm	Aesculap Hilan	3 rengas	Aesculap GB520R	Aesculap GB538R
NSXXXX.526	100.000 rpm	Aesculap Hilan XS	1 rengas	Aesculap GB740R	Aesculap GB741R
NSXXXX.525	100.000 rpm	Aesculap Hilan XS	2 rengas	Aesculap GB740R	Aesculap GB742R
NSXXXX.527	100.000 rpm	Aesculap Hilan XS	3 rengas	Aesculap GB740R	Aesculap GB743R
NSXXXX.515	100.000 rpm	Stryker	1 rengas	Stryker U2-DHil; 5100-100	Stryker „D-47“; 5100-10-47
NSXXXX.516	100.000 rpm	Stryker	2 rengas	Stryker U2-DHil; 5100-100	Stryker „D-48“; 5100-10-47

Figura numeru	Velocitate maxima	Sistem	Indicatore lungimea	Manipolo	Protezione dura madre
NSXXXX.025	20.000 giri/min	Aesculap Low Speed	N/D	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 giri/min	Aesculap Low Speed	N/D	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 giri/min	Aesculap Hilan	1 anello	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 giri/min	Aesculap Hilan	2 anelli	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 giri/min	Aesculap Hilan	3 anelli	Aesculap GB520R	Aesculap GB538R
NSXXXX.526	100.000 giri/min	Aesculap Hilan XS	1 anello	Aesculap GB740R	Aesculap GB741R
NSXXXX.525	100.000 giri/min	Aesculap Hilan XS	2 anelli	Aesculap GB740R	Aesculap GB742R
NSXXXX.527	100.000 giri/min	Aesculap Hilan XS	3 anelli	Aesculap GB740R	Aesculap GB743R
NSXXXX.515	100.000 giri/min	Stryker	1 anello	Stryker U2-DHil; 5100-100	Stryker „D-47“; 5100-10-47
NSXXXX.516	100.000 giri/min	Stryker	2 anelli	Stryker U2-DHil; 5100-100	Stryker „D-48“; 5100-10-47

Figurnummer	Maks. turtall	Sistem	Leigdøbetegnelse	Håndstykke	Dura mater-beskyttelse
NSXXXX.025	20.000 o/min	Aesculap Low Speed	N/A	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 o/min	Aesculap Low Speed	N/A	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 o/min	Aesculap Hilan	1-ring	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 o/min	Aesculap Hilan	2-ring	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 o/min	Aesculap Hilan	3-ring	Aesculap GB520R	Aesculap GB538R
NSXXXX.526	100.000 o/min	Aesculap Hilan XS	1-ring	Aesculap GB740R	Aesculap GB741R
NSXXXX.525	100.000 o/min	Aesculap Hilan XS	2-ring	Aesculap GB740R	Aesculap GB742R
NSXXXX.527	100.000 o/min	Aesculap Hilan XS	3-ring	Aesculap GB740R	Aesculap GB743R
NSXXXX.515	100.000 o/min	Stryker	1-ring	Stryker U2-DHil; 5100-100	Stryker „D-47“; 5100-10-47
NSXXXX.516	100.000 o/min	Stryker	2-ring	Stryker U2-DHil; 5100-100	Stryker „D-48“; 5100-10-47

Numărul figurii	Turație maximă	Sistem	Indicarea lungimii	Piesa de mână	Protecție dură
NSXXXX.025	20.000 U/min	Aesculap Low Speed	Nu se aplică	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 U/min	Aesculap Low Speed	Nu se aplică	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 U/min	Aesculap Hilan	1 inel	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 U/min	Aesculap Hilan	2 inele	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 U/min	Aesculap Hilan	3 inele	Aesculap GB520R	Aesculap GB538R
NSXXXX.526	100.000 U/min	Aesculap Hilan XS	1 inel	Aesculap GB740R	Aesculap GB741R
NSXXXX.525	100.000 U/min	Aesculap Hilan XS	2 inele	Aesculap GB740R	Aesculap GB742R
NSXXXX.527	100.000 U/min	Aesculap Hilan XS	3 inele	Aesculap GB740R	Aesculap GB743R
NSXXXX.515	100.000 U/min	Stryker	1 inel	Stryker U2-DHil; 5100-100	Stryker „D-47“; 5100-10-47
NSXXXX.516	100.000 U/min	Stryker	2 inele	Stryker U2-DHil; 5100-100	Stryker „D-48“; 5100-10-47

Številka	Maks. število vrtljajev	Sistem	Oznaka oblike	Ročnik	Zaščitna dura
NSXXXX.025	20.000 vrt./min	Aesculap Low Speed	/	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 vrt./min	Aesculap Low Speed	/	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 vrt./min	Aesculap Hilan	1 obroček	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 vrt./min	Aesculap Hilan	2 obročka	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 vrt./min	Aesculap Hilan	3 obročka	Aesculap GB520R	Aesculap GB538R
NSXXXX.526	100.000 vrt./min	Aesculap Hilan XS	1 obroček	Aesculap GB740R	Aesculap GB741R
NSXXXX.525	100.000 vrt./min	Aesculap Hilan XS	2 obročka	Aesculap GB740R	Aesculap GB742R
NSXXXX.527	100.000 vrt./min	Aesculap Hilan XS	3 obročka	Aesculap GB740R	Aesculap GB743R
NSXXXX.515	100.000 vrt./min	Stryker	1 obroček	Stryker U2-DHil; 5100-100	Stryker „D-47“; 5100-10-47
NSXXXX.516	100.000 vrt./min	Stryker	2 obročka	Stryker U2-DHil; 5100-100	Stryker „D-48“; 5100-10-47

Αριθμός συστήματος	Μεγ. αριθμός στροφών	Σύστημα	Χαρακτηριστικός δείκτης	Χειρολαβή	Προστασία της οπίσθιας μήτρας
NSXXXX.025	20.000 ο.α.λ.	Aesculap Low Speed	Δ/δωτ.	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 ο.α.λ.	Aesculap Low Speed	Δ/δωτ.	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 ο.α.λ.	Aesculap Hilan	1 δακτυλίου	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 ο.α.λ.	Aesculap Hilan	2 δακτυλίων	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 ο.α.λ.	Aesculap Hilan	3 δακτυλίων	Aesculap GB520R	Aesculap GB538R
NSXXXX.526	100.000 ο.α.λ.	Aesculap Hilan XS	1 δακτυλίου	Aesculap GB740R	Aesculap GB741R
NSXXXX.525	100.000 ο.α.λ.	Aesculap Hilan XS	2 δακτυλίων	Aesculap GB740R	Aesculap GB742R
NSXXXX.527	100.000 ο.α.λ.	Aesculap Hilan XS	3 δακτυλίων	Aesculap GB740R	Aesculap GB743R
NSXXXX.515	100.000 ο.α.λ.	Stryker	1 δακτυλίου	Stryker U2-DHil; 5100-100	Stryker „D-47“; 5100-10-47
NSXXXX.516	100.000 ο.α.λ.	Stryker	2 δακτυλίων	Stryker U2-DHil; 5100-100	Stryker „D-48“; 5100-10-47

Atēta numurs	Maks. apgrozāumu ātrums	Sistēma	Garuma apraksts	Uzgaļi	Smadzeņu cēlā aizsardzība
NSXXXX.025	20.000 U/min	Aesculap Low Speed	neatīstīcas	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 U/min	Aesculap Low Speed	neatīstīcas	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 U/min	Aesculap Hilan	1 gredzens	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 U/min	Aesculap Hilan	2 gredzeni	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 U/min	Aesculap Hilan	3 gredzeni	Aesculap GB520R	Aesculap GB538R
NSXXXX.526	100.000 U/min	Aesculap Hilan XS	1 gredzens	Aesculap GB740R	Aesculap GB741R
NSXXXX.525	100.000 U/min	Aesculap Hilan XS	2 gredzeni	Aesculap GB740R	Aesculap GB742R
NSXXXX.527	100.000 U/min	Aesculap Hilan XS	3 gredzeni	Aesculap GB740R	Aesculap GB743R
NSXXXX.515	100.000 U/min	Stryker	1 gredzens	Stryker U2-DHil; 5100-100	Stryker „D-47“; 5100-10-47
NSXXXX.516	100.000 U/min	Stryker	2 gredzeni	Stryker U2-DHil; 5100-100	Stryker „D-48“; 5100-10-47

Numer ury	Maks. predkioś obrotowa	System	Oznaczenie długości	Relejojśc	Ochrona opory twardej
NSXXXX.025	20.000 obr./min	Aesculap Low Speed	Nd.	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 obr./min	Aesculap Low Speed	Nd.	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 obr./min	Aesculap Hilan	1 pierścioń	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 obr./min	Aesculap Hilan	2 pierścioń	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 obr./min	Aesculap Hilan	3 pierścioń	Aesculap GB520R	Aesculap GB538R
NSXXXX.526	100.000 obr./min	Aesculap Hilan XS	1 pierścioń	Aesculap GB740R	Aesculap GB741R
NSXXXX.525	100.000 obr./min	Aesculap Hilan XS	2 pierścioń	Aesculap GB740R	Aesculap GB742R
NSXXXX.527	100.000 obr./min	Aesculap Hilan XS	3 pierścioń	Aesculap GB740R	Aesculap GB743R
NSXXXX.515	100.000 obr./min	Stryker	1 pierścioń	Stryker U2-DHil; 5100-100	Stryker „D-47“; 5100-10-47
NSXXXX.516	100.000 obr./min	Stryker	2 pierścioń	Stryker U2-DHil; 5100-100	Stryker „D-48“; 5100-10-47

Брój нутрашњеи	Макс. брój обротај	Систем	Означє дужине	Дршкє	Заштита тврє мозганє
NSXXXX.025	20.000 o/min	Aesculap Low Speed	није примењиво	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 o/min	Aesculap Low Speed	није примењиво	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 o/min	Aesculap Hilan	1 прстєн	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 o/min	Aesculap Hilan	2 прстєна	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 o/min	Aesculap Hilan	3 прстєна	Aesculap GB520R	Aesculap GB538R
NSXXXX.526	100.000 o/min	Aesculap Hilan XS	1 прстєн	Aesculap GB740R	Aesculap GB741R
NSXXXX.525	100.000 o/min	Aesculap Hilan XS	2 прстєна	Aesculap GB740R	Aesculap GB742R
NSXXXX.527	100.000 o/min	Aesculap Hilan XS	3 прстєна	Aesculap GB740R	Aesculap GB743R
NSXXXX.515	100.000 o/min	Stryker	1 прстєн	Stryker U2-DHil; 5100-100	Stryker „D-47“; 5100-10-47
NSXXXX.516	100.000 o/min	Stryker	2 прстєна	Stryker U2-DHil; 5100-100	Stryker „D-48“; 5100-10-47

Figur nummer	Max. varvantal	System	Långbeteckningar	Handstycke	Skydd av dura mater
NSXXXX.025	20.000 V/min	Aesculap Low Speed	n/a	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 V/min	Aesculap Low Speed	n/a	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 V/min	Aesculap Hilan	1 ring	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 V/min	Aesculap Hilan	2 ringar	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 V/min	Aesculap Hilan	3 ringar	Aesculap GB520R	Aesculap GB538R
NSXXXX.526	100.000 V/min	Aesculap Hilan XS	1 ring	Aesculap GB740R	Aesculap GB741R
NSXXXX.525	100.000 V/min	Aesculap Hilan XS	2 ringar	Aesculap GB740R	Aesculap GB742R
NSXXXX.527	100.000 V/min	Aesculap Hilan XS	3 ringar	Aesculap GB740R	Aesculap GB743R
NSXXXX.515	100.000 V/min	Stryker	1 ring	Stryker U2-DHil; 5100-100	Stryker „D-47“; 5100-10-47
NSXXXX.516	100.000 V/min	Stryker	2 ringar	Stryker U2-DHil; 5100-100	Stryker „D-48“; 5100-10-47

Ábrázám	Max. fordulatszám	Rendszer	Hosszúzájelölés	Kézidrab	Duravédelem
NSXXXX.025	20.000 ford./perc	Aesculap Low Speed	n.a.	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 ford./perc	Aesculap Low Speed	n.a.	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 ford./perc	Aesculap Hilan	1 gyűrűs	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 ford./perc	Aesculap Hilan	2 gyűrűs	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 ford./perc	Aesculap Hilan	3 gyűrűs	Aesculap GB520R	Aesculap GB538R
NSXXXX.526	100.000 ford./perc	Aesculap Hilan XS	1 gyűrűs	Aesculap GB740R	Aesculap GB741R
NSXXXX.525	100.000 ford./perc	Aesculap Hilan XS	2 gyűrűs	Aesculap GB740R	Aesculap GB742R
NSXXXX.527	100.000 ford./perc	Aesculap Hilan XS	3 gyűrűs	Aesculap GB740R	Aesculap GB743R
NSXXXX.515	100.000 ford./perc	Stryker	1 gyűrűs	Stryker U2-DHil; 5100-100	Stryker „D-47“; 5100-10-47
NSXXXX.516	100.000 ford./perc	Stryker	2 gyűrűs	Stryker U2-DHil; 5100-100	Stryker „D-48“; 5100-10-47

Paveikėlijo numeris	Maks. apskiuo skaičius	Sistema	Ilgio pavidėnimas	Rankinio falais	Kietjo dangalo apsauga
NSXXXX.025	20.000 aps./min	Aesculap Low Speed	N/A	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 aps./min	Aesculap Low Speed	N/A	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 aps./min	Aesculap Hilan	1 žiedo	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 aps./min	Aesculap Hilan	2 žedų	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 aps./min	Aesculap Hilan	3 Ring	Aesculap GB520R	Aesculap GB538R
NSXXXX.526	100.000 aps./min	Aesculap Hilan XS	1 žiedo	Aesculap GB740R	Aesculap GB741R
NSXXXX.525	100.000 aps./min	Aesculap Hilan XS	2 žedų	Aesculap GB740R	Aesculap GB742R
NSXXXX.527	100.000 aps./min	Aesculap Hilan XS	3 Ring	Aesculap GB740R	Aesculap GB743R
NSXXXX.515	100.000 aps./min	Stryker	1 žiedo	Stryker U2-DHil; 5100-100	Stryker „D-47“; 5100-10-47
NSXXXX.516	100.000 aps./min	Stryker	2 žedų	Stryker U2-DHil; 5100-100	Stryker „D-48“; 5100-10-47

Número de figura	Rotações máximas	Sistema	Designação do comprimento	Peça de mão	Proteção dural
NSXXXX.025	20.000 rpm	Aesculap Low Speed	N/A	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 rpm	Aesculap Low Speed	N/A	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 rpm	Aesculap Hilan	1 anel	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 rpm	Aesculap Hilan	2 anéis	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 rpm	Aesculap Hilan	3 anéis	Aesculap GB520R	Aesculap GB538R
NSXXXX.526	100.000 rpm	Aesculap Hilan XS	1 anel	Aesculap GB740R	Aesculap GB741R
NSXXXX.525	100.000 rpm	Aesculap Hilan XS	2 anéis	Aesculap GB740R	Aesculap GB742R
NSXXXX.527	100.000 rpm	Aesculap Hilan XS	3 anéis	Aesculap GB740R	Aesculap GB743R
NSXXXX.515	100.000 rpm	Stryker	1 anel	Stryker U2-DHil; 5100-100	Stryker „D-47“; 5100-10-47
NSXXXX.516	100.000 rpm	Stryker	2 anéis	Stryker U2-DHil; 5100-100	Stryker „D-48“; 5100-10-47

Figurálne číslo	Max. otáčky	System	Označení délky	Rukovět	Ochrana pred epidurálnym hematómou
NSXXXX.025	20.000 ot./min	Aesculap Low Speed	Nie je k dispozíci	Aesculap GB267R	Aesculap GB266M
NSXXXX.25A	20.000 ot./min	Aesculap Low Speed	Nie je k dispozíci	Aesculap GB267R	Aesculap GB266M
NSXXXX.512	100.000 ot./min	Aesculap Hilan	jednokruhové	Aesculap GB520R	Aesculap GB536R
NSXXXX.513	100.000 ot./min	Aesculap Hilan	dvojkruhové	Aesculap GB520R	Aesculap GB537R
NSXXXX.514	100.000 ot./min	Aesculap Hilan	trojkruhové	Aesculap GB520R	Aesculap GB538R